

TABLE OF CONTENTS

Mechanical and Electronic Counters 2

Mechanical and Electronic Counters



Counters, manual four figures, with button for manual counting

Item number:
40344

The casing made of chrome-plated steel lies comfortably in the hand. The manual counter is being held with the forefinger at the finger bow. The counter is turned back to zero by using the lateral reset button.

Item No.	Counting limit
40344010	0 - 9999

[» go to product \(website\)](#)



Blood cell counters

Item number:
40345

Electronic counters and memory devices for microscopic blood testing

The electronic memory devices specify blood counting. The microscope image can be observed permanently, the numbers can be typed in easily. Every device comes with an instruction manual.

Advantages etc.:

- Power supply operation
- Rechargeable battery operation
- Software (download) with computer analysis
- Big TFT color display
- Storage of last counting on device

Interfaces

Data transfer - USB mini

Power supply

Power supply operation - USB adapter 5VDC / 300mA

Rechargeable - USB power bank external
battery operation

Functional keys

Labelling On/Off, -/Eval, +, Go On, Corr. Clear

Counting keys

Labelling

40345010 - Eos, Baso, Promy, Myelo, Div, Juve, Stab, Seg, Lym, Mono

40345015 - At Ly, My Bl, Met My, Myelo, Normo, Lym, Seg, Mono, Plas, Tox Gr, Stab, Eos, Baso, Div/Reti, Div/Ery

Programs

Counting steps - 10-1000 in 10 steps

Programs - LEUKO, MYELO/RBC, RETI, FREE

Item No.	Number of counting keys	Dimensions	Number of functional keys	Colour	Counting limit	Display
40345010	10	130 x 235 x 35 mm	5	blue	max. 1000	TFT color display 3,5'
40345015	15	130 x 235 x 35 mm	5	blue	max. 1000	TFT color display 3,5'

[» go to product \(website\)](#)

Index

B

Blood cell counters 2

C

Counters, manual four figures, with button for manual counting 1